



# LIP PLUMPER

USER GUIDE





## BE CONFIDENT IN YOUR COSMETICS

L(A)B is the creation of Oxford University Scientists and Sisters Lauren, Charlotte and Sophie with the aim to improve the effectiveness of your everyday skincare and healthcare regimes through cutting edge scientific devices that fit easily into your life.

L(A)B is committed to delivering products of the highest innovation and quality from concept to packaging. We are a unique brand, the first British, Family run Electronic Skincare + Healthcare company. L(A)B are the Beauty Tool Experts.

## For fuller lips

5 technologies in one device specifically designed to target thinning lips, loss of fullness and the lines and wrinkles around the lips.

- Radio Frequency
- Micro-current
- Galvanic Current
- Thermal Therapy
- Massage



# How to Use

1. Remove any lip gloss or lip stick
2. Apply the lip volumising balm to the lips, and the edges of mouth over any lines and wrinkles.
3. Turn on the lip plumper by pressing and holding the button
4. Place the treatment head to the lips and treat all the lips and also the lines around the lips on to the skin using the following 2 motions for **2 minutes twice a day**:
  - Move the lip plumper over the lips like you are applying a lipstick. Start in the centre of the lips and move to the edges in steady motions. Ensure you cover the whole lip area, and also the skin around the lips where the fines and wrinkles are
  - Use small circular motions across the lips and the skin around the lips. Start in the centre and work outwards.

## How to Charge

- 1.** Insert the charging cable in to the charging point on the device, and connect the other end to a USB charging point.
- 2.** The indicator light flashes blue when charging is required and is a constant blue light when fully charged.
- 3.** A full charge takes 2 hours and will last up to 12 hours.

## Cautions

- 1.** The device does include a heating function and will feel warm.
- 2.** If you feel any discomfort stop use immediately.
- 3.** The device is not waterproof, do not submerge in water.
- 4.** Store in a dry well ventilated place and avoid direct sun light.
- 5.** Do not use if you have an electrical implant fitted such as a pace maker.
- 6.** Do not use if you are pregnant.





Contact us

E: [info@lab-lifeandbeauty.com](mailto:info@lab-lifeandbeauty.com) | [www.lab-lifeandbeauty.com](http://www.lab-lifeandbeauty.com)

