

CelluFade BodyBliss

324-317

How to Use:

Body Mode:

While in "BODY" setting, both the LIPO Mode (Vacuum Therapy) and the EMS Mode (Electrical Muscle Stimulation) can be used either together or separately with 3 different intensity levels for each.

While using LIPO Mode - Ensure to first apply a gel or cream to the treatment area and keep your finger on the RELEASE button in the center. Do not use the LIPO feature directly on one area for longer than 5 seconds.

Attachments:

Use the big suction cup for larger areas such as legs, stomach, and buttocks.

Use the small suction cup for smaller areas such as arms.

Face Mode:

While in "FACE" setting, only the EMS Mode (Electrical Muscle Stimulation) is available for use.

LIPO mode cannot be used on the face; no attachment should be used for when performing a face treatment.

Adjust the intensity level per your comfort and desired results.

The EMS function is not to be used directly under the eye area.

Make sure to apply a cream, gel or serum on the treatment area to allow for optimal conductivity of the EMS technology.