

LaserBoost Advanced RF Hair Growth Brush

324-320

How to Use:

Cleanse hair and scalp prior to use, hair may be damp or dry during treatment.

After turning on the device by tapping the power button twice, select desired mode (RF / Laser / Vibration) and intensity levels.

- RF Radio Frequency Mode: 5 Intensity Levels
- Laser Mode: No Intensity Level – this is a light-emitting mode
- Vibration Massaging Mode: 5 Intensity Levels

Gently glide the LaserBoost+ from the front to the back of your scalp in slow combing motions. Repeat for the duration of the 5-minute treatment.

The device will emit a sound every couple of seconds (length varies with each mode) as a reminder to move the device to a different area.

The device will automatically shut off after 15 minutes.

For best results, it is recommended to combine all 3 technologies during treatment.