



SUPREME SMILE

TEETH WHITENING KIT

USER MANUAL



ABOUT

The Supreme Smile Teeth Whitening Kit is the ideal solution for whitening any discoloration resulting from: aging, coffee, tea, red wine, tobacco or food coloring. The device utilizes cold LED light technology combined with whitening gel for an outstanding result. This advanced whitening process is delivered through the mouthpiece. Applying the principles of dental engineering, the soft mouthpiece offers optimal comfort for a truly seamless whitening process. Prior to use, please read this instruction manual thoroughly.

LIGHT MODES:

- **Blue Mode-** the rays emitted from the cold blue LED light react with the teeth whitening gel to remove stains.
- **Red Mode-** the rays emitted from the red LED light eliminate bacteria and simultaneously soothe your gums. Red light irradiation is believed to strengthen cell metabolism, which can help promote the healing of wounds and ulcers. It also helps increase the phagocytosis of white blood cells, improving your teeth's immune function.

INCLUDES



LED Light



Results Charts



Pens



Instruction Manual

WHITENING STEPS

1. Rinse both your mouth and whitening mouthpiece.
2. Compare the color of your teeth against the provided shade guide and record your stage.
3. Twist whitening pen until gel begins to flow. Apply gel evenly to the lower and upper teeth. Please note, do not apply too much gel, it will not speed up the process- each use should only be: .5-.7ml.
4. Insert mouthpiece and bite down.
5. Turn on the device light.
6. To avoid diluting the gel, take a deep breath at the start of treatment for minimized saliva production.
7. The device will automatically switch off after 15 minutes on blue mode and 10 minutes on blue/red mode. When treatment finishes, thoroughly rinse both your mouth and the mouthpiece.
8. Compare your teeth again using the shade guide to assess your progress.
9. Do not eat or drink anything except water for 30 minutes following treatment.

WHITENING STEPS: QUICK GUIDE



Rinse mouth with water



Apply the gel on tooth surface



Put the mouthpiece into the mouth and bite it, then turn on the light



15 minute treatment (Blue mode)
10 minute treatment (Blue&Red mode)



Rinse mouth with water



Rinse mouthpiece

PRECAUTIONS AND ADDITIONAL INFORMATION

1. Avoid brushing gel on gums or the inner mouth/lips. While not serious, contact will result in a burning feeling. If this happens, rinse mouth with water and the sensation will disappear within 24 hours.
2. Do not use if allergic to gel.
3. Do not swallow the gel.
4. Keep out of reach of children.
5. Do not use if pregnant or breastfeeding.
6. Should the gel get on clothing, wash immediately.
7. If the gel gets on your face or skin, wash immediately.
8. To preserve the gel for future uses, always replace cap on gel pen.
9. Store in a cool, dry area avoiding sunlight.
10. Not for use on anyone under 18.
11. Do not use if you have sensitive teeth.
12. Do not use on any decaying or loose teeth.

** Please Note: This device will not whiten any existing restoration such as: fillings, veneers, bridges, or crowns.

INGREDIENTS

Carbamide Peroxide, Glycerol, Propylene Glycol, Deionized Water, Carbomer, Carboxymethyl, Poly Vinypyrrolidone, Sodium Hydroxide, Menthol.

STORAGE

Avoid sunlight and keep in cool and dry area. The gel pens can and should be kept in refrigerator.

SHELF LIFE

12 months

LIMITED LIFETIME WARRANTY

Your Supreme Smile Teeth Whitening Kit is backed by a limited lifetime manufacturer's warranty. Supreme Smile will repair or replace your device at any time should it fail due to a defect in material or workmanship, subject to the certain limitations.

This limited warranty does not cover any damage that results from unauthorized or improper use, service, or repair. Further, it does not cover damage caused by accident, impact, negligence, or normal wear and tear. Should you discover your Supreme Smile Teeth Whitening Kit is not

functioning properly, please send your device to our repair center for evaluation. If your product cannot be repaired or serviced, we will reserve the right to change it for a similar or newer model.

Please note that a flat rate of \$35.00 will be charged to cover a service evaluation fee and return shipping of your device. All warranty claims must be accompanied by a copy of your proof of purchase from an authorized retailer. Please send your device, proof of purchase, and a check or money order in the amount of \$35.00 made out to Supreme Smile to:

Address:

Supreme Smile Service Center
3069 Taft Street
Hollywood, FL 33021

Contact:

warranty@supreme-smile.com